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*The Fuego SCS basic represents a new standard among laboratory gas burners. With its innovative design this laboratory gas burner is suitable for all flame related applications in the laboratory.*

## Modern Design

*The unique shape is eye-catching\*. The low profile body facilitates ergonomic operation; the streamlined design reduces air flow disturbances in a clean room workbench to a minimum. Moreover, the Fuego SCS basic is very easy to clean – a laboratory gas burner without grooves and rough edges. Clogging is not an issue. Spilled liquids simply escape through an integrated drain while solid substances can easily be removed by opening the burner shaft. It is even possible to dismantle the individual components of the removable burner head for in-depth cleaning.*

## Top Rate Safety

*The Safety Control System SCS® demonstrates state-of-the-art safety technology which constantly analyzes potential hazards and, if necessary, initiates safety measures, such as an interruption of the gas supply. In addition to the ignition and flame control function and the overheating protection, the new safety concept also features continuous burner head control (BHC). BHC unfailingly identifies burner head clogging by liquids or solid*

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Fuego SCS basic  
Safety laboratory gas burner  
with button function and foot pedal

Art.-Nr./Art.-No.  
8.001.000

*substances and ensures the correct assembly of the burner head. Exceptional passive safety features: After prolonged use, the residual heat display will signal that the burner head is still hot in order to protect the user from burns.*

## Simple and Safe Operation

*The Fuego SCS basic can be operated with a foot pedal or button function. Different programs are available: a flexible start-stop function and the conventional foot pedal operation for short flame sterilization. Additionally the flame can be started and stopped easily by pushing the function knob.*

## Robust

*Fabricated entirely of stainless steel, with fireproof controls and a display, protected by heat resistant glass, the Fuego SCS basic can withstand extreme lab conditions.*

**Fuego SCS basic - Ideal for use in cleanroom workbenches and laboratories.**

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## The range:

### Fuego SCS basic Art.-No. 8.001.000

with button function and stainless steel foot pedal  
3 standard-programs for button (function knob) and foot pedal  
SCS (Safety Control System) with BHC (Burner Head Control)  
Removable and decomposable burner head  
Tilt mechanism, right / left  
Holding device for 3 inoculation loop holders  
Nozzles for natural gas, propane / butane gas  
Turbo flame  
Wrench 17 mm for gas connection  
Screwdriver for burner head and cover of burner shaft  
Tubing connector with swivel nut  
Switching power supply (global)  
Instruction manual and 2-year warranty

#### Technical Data:

<b>Technology</b>	Microprocessor
<b>Programs</b>	
Foot pedal:	Standard (flame during pressed foot pedal) Start-Stop with timer, 60 min Start-Stop with timer, 60 min
Button:	Start-Stop with timer, 60 min

#### Safety features

Safety Control System (SCS)  
with gas safety cut off: ignition and flame control  
temperature monitor  
burner head clogging  
and assembly monitor (BHC)  
automatic unit switch off, 4 h  
residual heat display

#### Gas supply and consumption

Gas connection: 1/4" left with gas filter  
Gas types: I12ELL3B/P: natural gas E/LL, 18 - 25 mbar  
liquid gas, 20 - 50 mbar  
Connected load: 70 g/h liquid gas  
Continuous cartridge operation: CV 360 - 40 min, Express 444 - 50 min,  
CG 1750 - 150 min, C 206 - 170 min,  
CP 250 - 210 min, CV 470 - 370 min

#### Temperatures

Flame temperature: 1350 °C on liquid gas  
1300 °C on natural gas (E)  
Temperature threshold level: 1.0 kW liquid gas, 1.0 kW natural gas

#### Electrical


Power consumption: 2 VA  
Power connection: 100 - 240 V / 50/60 Hz / max. 0.3 A  
9 V DC / 1 A

#### Mechanical

Casing and operating controls: stainless steel / glass, UV and solvent resistant  
Burner head: removable and decomposable, stainless steel  
Cover of burner shaft: Ø 23 mm, with drains  
Measurements (w x h x d): 103 x 49 x 130 mm  
Weight: 700 g

#### Licenses

DIN-DVGW Reg.-No.: NG-2211AS0167  
CE: EN 61326-1, EN61000-3-2, EN 61010  
EEC guidelines: 89/336/EEC, 73/23/EEC

**Turboflamme**  
im Handumdrehen eine extra  
harte und stabile Flamme  
**Turbo flame**  
an extremely firm and  
consistent flame with only a few  
flicks of the wrist

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**Fußpedal**  
aus Edelstahl  
**Foot pedal**  
made of  
stainless steel

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**WLD-TEC** GmbH