

## Thermo Scientific Fiberlite F9-4x1000y Rotor

A blue Thermo Scientific F8-4x1000 centrifuge rotor. The rotor is cylindrical with a blue body and a silver-colored top rim. The top surface is black with four large, curved, black slots for tubes. The Thermo Scientific logo and model number are printed on the front.

<b>Type</b>	Fixed Angle
<b>Material</b>	Carbon Fiber Composite
<b>Capacity</b>	4 x 1000 mL
<b>Maximum Speed</b>	9,000 rpm
<b>Tube Dimensions Ø x L</b>	98 x 95 mm
<b>Angle</b>	22°
<b>K Factor</b>	5,950
<b>Net Weight</b>	14.29 kg/31.5 lbs
<b>Warranty</b>	15 years

	RCF	RADIUS
Maximum	15,241 xg	168 mm
Minimum	2,268 xg	25 mm

		THERMO SCIENTIFIC			
		SORVALL RC-6 Plus/Evolution RC Series	SORVALL RC-6/RC-5/ RC-26/RC-2	JOUAN KR22i	JOUAN KR25i
Fiberlite F9-4x1000y	Cat. No.	76981	76981	096-041042	096-041041
	Max. Speed	8,000/9,000 rpm	7,000 rpm	9,000 rpm	9,000 rpm
	Max. RCF	12,042/15,241 xg	9,220 xg	15,241 xg	15,241 xg

Rotor Lid Assembly	099-041004	099-041004	099-041042	099-041041
O-ring Refresher Kit <sup>1</sup>	021-041004	021-041004	021-041042	021-041041

									REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	Type	Cat. No.	Qty. /Pkg.	Description
1000.0	1000.0	Nalgene® PPCO Oak Ridge Bottle	3141-1002	2	4	9,000	15,241	98 x 195	Closure	Included	2	PP Sealing
1000.0	1000.0	Nalgene PC Oak Ridge Bottle	3140-1002	2	4	9,000	15,241	98 x 195	Closure	Included	2	PP Sealing
500.0	450.0	Nalgene PPCO Oak Ridge Bottle	3141-0500	24	4	9,000	15,241	70 x 160	Closure Adapter	Included 010-0145	24 2	PP Sealing 1 place/adapter
500.0	450.0	Nalgene PC Oak Ridge Bottle	3140-0500	24	4	9,000	15,241	70 x 160	Closure Adapter	Included 010-0145	24 2	PP Sealing 1 place/adapter
250.0	—	Corning® Disposable Conical Tube²	—	—	4	9,000	—	—	Adapter	010-1096	2	1 place/adapter
250.0	250.0	Nalgene PPCO Oak Ridge Bottle	3141-0250	36	4	9,000	15,241	61 x 133	Closure Adapter	Included 010-0150	— 2	PP Sealing 1 place/adapter
250.0	250.0	Nalgene PC Oak Ridge Bottle	3140-0250	36	4	9,000	15,241	61 x 135	Closure Adapter	Included 010-0150	36 2	PP Sealing 1 place/adapter

**PERFECT FIT: FREE  
SET OF NALGENE  
PRODUCTS  
INCLUDED WITH  
ROTOR PURCHASE**

<sup>2</sup> Performance depends upon tube brand. Contact your sales representative for recommended tube manufacturers and further information.

## LABWARE ORDERING INFORMATION CONTINUED

Tube Vol. (mL)	Vol. Fill (mL)	Description	Cat. No.	Qty. /Pkg.	Tubes /Rotor	Max. Speed (rpm)	Max. RCF (xg)	Dims ØxL (mm)	REQUIRED TUBE CLOSURES, ADAPTERS, TOOLS & ACCESSORIES			
									Type	Cat. No.	Qty. /Pkg.	Description
175.0	175.0	Nalgene PPCO Conical Wide-Mouth Tube	3143-0175	36	4	9,000	15,241	61 x 144	Closure Adapter	Included 010-1132	36 2	PP Sealing 1 place/adapter
175.0	175.0	Nalgene PC Conical Wide-Mouth Tube	3144-0175	36	4	9,000	15,241	61 x 144	Closure Adapter	Included 010-1132	36 2	PP Sealing 1 place/adapter
80.0	73.0	PC Oak Ridge Tube	010-0515	6	12	9,000	15,241	38 x 109	Closure Adapter	Included 010-1093	6 2	PP Sealing 3 place/adapter
80.0	76.0	PPCO Oak Ridge Tube	010-1280	6	12	9,000	15,241	38 x 109	Closure Adapter	Included 010-1093	6 2	PP Sealing 3 place/adapter
85.0	81.0	Nalgene PC Oak Ridge Tube	3118-0085	100	12	9,000	15,241	38 x 109	Closure Adapter	Included 010-1093	100 2	PP Screw Top 3 place/adapter
50.0	50.0	Nunc™ Disposable Conical Tube <sup>2</sup>	334959	25	20	9,000	15,241	30 x 115	Closure Adapter	Included 010-0180	25 2	PP Sealing 5 places/adapter
50.0	—	Corning Disposable Conical Tube <sup>2</sup>	—	—	20	9,000	15,241	—	Adapter	010-0180	2	5 places/adapter
50.0	—	BD Falcon® Disposable Conical Tube <sup>2</sup>	—	—	20	9,000	15,241	—	Adapter	010-0180	2	5 places/adapter
50.0	46.0	Nalgene PC Oak Ridge Tube	3138-0050	50	24	9,000	15,241	29 x 115	Closure Adapter	Included 010-1091	50 2	PP Sealing 6 places/adapter
50.0	46.0	Nalgene PPCO Oak Ridge Tube	3139-0050	50	24	9,000	15,241	29 x 115	Closure Adapter	Included 010-1091	50 2	PP Sealing 6 places/adapter
50.0	46.0	Nalgene FEP Oak Ridge Tube	3114-0050	10	24	9,000	15,241	29 x 108	Closure Adapter	Included 010-1091	10 2	ETFE Screw Cap 6 places/adapter
30.0	30.0	Nalgene PC Oak Ridge Tube	3138-0030	50	28	9,000	15,241	26 x 102	Closure Adapter	Incl 010-		
30.0	30.0	Nalgene PPCO Oak Ridge Tube	3139-0030	50	28	9,000	15,241	26 x 102	Closure Adapter	Incl 010-		
30.0	30.0	Nalgene FEP Oak Ridge Tube	3114-0030	10	28	9,000	15,241	26 x 102	Closure Adapter	Incl 010-		
16.0	16.0	Nalgene PC Oak Ridge Tube	3139-0016	50	56	9,000	15,241	18 x 107	Closure Adapter	Incl 010-		
16.0	16.0	Nalgene PPCO Oak Ridge Tube	3138-0016	50	56	9,000	15,241	18 x 107	Closure Adapter	Incl 010-		
15.0	—	Corning Disposable Conical Tube <sup>2</sup>	—	—	48	—	—	—	Adapter	010-		
15.0	—	BD Falcon Disposable Conical Tube <sup>2</sup>	—	—	48	—	—	—	Adapter	010-		
10.0	8.0	Nalgene PPCO Oak Ridge Tube	3139-0010	50	72	9,000	15,241	16 x 82	Closure Adapter	Included 010-1307	50 2	PP Sealing 18 places/adapter
10.0	8.0	Nalgene PC Oak Ridge Tube	3138-0010	50	72	9,000	15,241	16 x 82	Closure Adapter	Included 010-1307	50 2	PP Sealing 18 places/adapter
10.0	—	Vacutainer® Tube <sup>2</sup>	—	—	72	—	—	16 x 100	Adapter	010-1085	2	18 places/adapter
3.0	—	Vacutainer Tube <sup>2</sup>	—	—	120	—	—	10 x 64	Adapter	010-1131	2	30 places/adapter

<sup>1</sup> Includes 5 O-rings, 3 O-ring Grease Packets and 4 oz FiberClene Lab Cleaner (017-0024).

<sup>2</sup> Performance depends upon tube brand. Contact your sales representative for recommended tube manufacturers and further information.

Note: Thermo Scientific Sorvall RC-6 Plus centrifuges manufactured prior to March 2009 may require software upgrade 314295G01 (NA) or HKKS308102E (Europe and Asia Pacific).

© 2010 Thermo Fisher Scientific Inc. All rights reserved. Beckman is a registered trademark of Beckman Coulter Inc. Hitachi is a registered trademark of Nissei Sangyo America. Corning is a registered trademark of Corning Inc. BD Falcon and Vacutainer are registered trademarks of BD Bioscience. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For more information, please visit  
[www.thermoscientific.com/rotors](http://www.thermoscientific.com/rotors).

**Thermo**  
 SCIENTIFIC